

Collagen Solution for Tissue Engineering Protocol

I. Compatible Products

Product

 UltraPure collagen I for tissue engineering, bovine placenta (001-G70-MF1)

 UltraPure collagen III for tissue engineering, bovine placenta (001-G70-MF3)

 Collagen I for tissue engineering, bovine placenta (001-G70-MF5)

 Collagen I for tissue engineering, bovine skin (001-G70-MF7)

 UltraPure collagen I for tissue engineering, human placenta (009-G70-MF1)

I. Materials

Product

 Collagen: protein concentration is around 2 mg/mL

 5X cell culture medium (RPMI, M199 or DMEM)

 10X PBS (sterile filtered)

 Neutralization Solution (0.1 M NaOH, sterile filtered)

II. Procedure for Collagen Gel Solution Preparation

Note: Keep all solutions on ice during the solution preparation.

1. According to Table 1, prepare the desired volume of collagen gel solution in the following sequence.

Table 1. Preparation of collagen gel solution

Total Volume (mL)	Collagen Solution (mL)	Medium (mL)	Neutralization Solution (μL)
10	8	2	250
5	4	1	125
2.5	2	0.5	62.5
1	0.8	0.2	25
0.75	0.6	0.15	18.75
0.5	0.4	0.1	12.5

Note: 5X medium should be compatible with the desired cell type.

a.) In a sterile tube, add the appropriate volume of collagen solution.

b.) Next, add the corresponding volume of 5X PBS or medium. Mix well.

Note: Solution should be yellow in color

c.) Finally, add the Neutralization Solution and immediately mix well.

Note: Solution should change to pink/red color

2. After mixing, keep the solution on ice. The pH of the collagen gel solution should be neutral, which is indicated by the pink/red color of phenol red in the 5X medium.

III. Procedure for Top Culture

1. Harvest and resuspend desired cell line at 0.1 to 2.0 x 10⁶ cells/mL.
2. Pipette a proper size of a chilled collagen gel solution onto a tissue culture plate or dish. See Table 2 for recommended measurements.
3. Immediately transfer to 37°C incubator for 60 minutes to initiate polymerization of the collagen.

Note: The polymerized gel will look cloudy.

Table 2. Suggested collagen gel amounts

Culture Dish	Volume (mL)
96-well	0.1
48-well	0.2
24-well	0.5
6-well	1
3-mm	2
60-mm	3
100-mm	5

4. After formation of the collagen gel, seed desired cells onto the collagen gel.
5. Overlay polymerized collagen gel with culture media.
6. Incubate cells overnight or several days at 37°C with appropriate CO₂ levels. Change medium daily.
7. Cells can be visualized using phase contrast microscopy and can be analyzed by other tests such as cell viability assay.

Rockland products are for research use only and are not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Immunochemicals, Inc., P.O. BOX 5199, Limerick, Pennsylvania 19468, USA.

© 2018 Rockland Immunochemicals, Inc. All rights reserved.