

Protocol

DYKDDDDK (FLAG® Tag) Immunoprecipitation

Anti-DYKDDDDK Affinity Gel is a purified mouse IgG2a monoclonal antibody coupled to activated agarose. This product is intended for purification of proteins containing the FLAG® epitope tag sequence.

Anti-DYKDDDDK Affinity Gel is optimally suited for immunoprecipitation and purification of FLAG® tagged proteins. Anti-DYKDDDDK Affinity Gel antibody recognizes the FLAG® epitope tag fused to the amino terminal, carboxyl terminal ends or an internal location of targeted fusion proteins. The epitope tag peptide sequence was first derived from the 11-amino-acid leader peptide of the gene-10 product from bacteriophage T7. DYKDDDDK is the most commonly used hydrophilic octapeptide tag. The control peptide D-Y-K-D-D-D-K is recommended for competitive elution to recover proteins.

Note: Remove only the amount of gel necessary for the intended purification (binding capacity is >1 µg protein per 10 µL of gel). Swirl the vial to thoroughly resuspend the gel, remove the desired amount, and transfer the gel to a microcentrifuge tube.

I. Reagents Required

Product	Item No.
DYKDDDDK (FLAG®) Peptide	000-000-383
Wash Buffer: 10X TBS pH 7.5	MB-012-0100
DYKDDDDK (FLAG®) Positive Control Lysate	W00-001-383
Triple FLAG® Peptide	000-000-H47
2X SDS-PAGE Sample Buffer without DTT or b-ME	MB-018
DYKDDDDK Tag (Anti-FLAG®) Antibody Peroxidase Conjugated	200-303-383

Note: Anti-DYKDDDDK Affinity Gel is supplied in 0.02M Potassium Phosphate, 0.15M Sodium Chloride pH 7.2 in 50% (v/v) glycerol and 0.01% (w/v) sodium azide (added as a preservative). Store vial at 4° C prior to opening vial. The unopened product is stable for one year when stored appropriately. After use, anti-DYKDDDDK Affinity Gel should be cleaned and stored in with PBS and azide. Do not freeze in the absence of glycerol.

II. Binding Specificity

Anti-DYKDDDDK Affinity Gel binds the FLAG® epitope tag sequence (Asp-Tyr-Lys-Asp-Asp-Asp-Lys) fused to the amino terminal, carboxy terminal or internal locations of targeted recombinant proteins expressed in transfected or transformed cells. D-Y-K-D-D-D-K peptide is recommended for competitive elution to recover fusion protein.

III. Procedure

1. Clarify the extract containing the FLAG® tagged protein obtained from cell lysates (provided by user) to remove any insoluble material.
Note: Typically insoluble material is removed by centrifugation at 12,000 rpm for 15-20 minutes, or other similar methods.
2. Filter the protein extract using a 0.22 µm or 0.45 µm filter to remove any remaining particulate that otherwise may obstruct the flow of buffer through the gel in later steps.
3. Transfer the protein extract to the microcentrifuge tube containing the Anti-DYKDDDDK Affinity Gel and incubate for at least 3 hours at 4°C with gentle mixing to capture the FLAG® tagged proteins.
Note: Mixing should be done using a rocking platform. Do not use a magnetic stirring system because this will destroy the gel beads. Some proteins may require overnight incubations to maximize the opportunity for antibody to bind. In these instances, it may be necessary to add protease inhibitors and/or bacteriostatic agents (provided by user) to prevent proteolysis and microbial growth.
4. Once the binding step is complete, collect the affinity gel from the microcentrifuge tube by centrifugation (2000-3000 rpm for 5 minutes) or by filtration using an appropriately sized apparatus.
5. Wash the affinity gel with 10 gel volumes of Wash Buffer to remove binding of nonspecific proteins. Check the last wash by UV absorbance to confirm that no protein is present. Repeat washing if needed.
6. After the last centrifugation, aspirate the Wash Buffer from the affinity gel.
7. Elute FLAG® tagged protein from the affinity gel by peptide competition using the control FLAG® peptide D-Y-K-D-D-D-K (purchased separately).
8. Add 3 gel volumes of control peptide to the affinity gel prepared from a stock concentration to a final concentration of 0.5 mg/mL in Wash Buffer.
9. Incubate the mixture for 30 minutes at 4°C, mixing as before.
10. Centrifuge to collect the affinity gel and recover the supernatant, which contains the eluted FLAG® tagged protein.
Note: The amount of FLAG® tagged protein present can be monitored by measuring the absorbance of the supernatant at 280 nm (if the extinction coefficient for the protein is known). Elution conditions using Triple FLAG® peptide must be optimized by the user.

Notes

- Larger-scale purification may be optimized for column chromatography. In these instances, protease inhibitors and/or bacteriostatic agents should be added to Wash Buffer.
- Conditions for using SDS-PAGE Sample Buffer, 0.1 M Glycine, pH 3.5, or other high pH or low pH elution buffers to elute FLAG® tagged protein from the affinity gel must be optimized by the user. Exposure to harsh conditions may inactivate the ability of the antibody to bind protein thereby inactivating the gel for reuse.
- Do not use the following buffers:
 - a. Chaotropic reagents such as guanidine HCl
 - b. Urea or sodium chloride greater than 1M
 - c. Reducing agents DTT, DTE, or 2-mercaptoethanol
 - d. TWEEN 20 or TRITON X-100 greater than 5%
 - e. IGEPAL CA-630 or CHAPS greater than 0.1%

- f. Digitonin greater than 0.2%
 - g. Deoxycholate or SDS
11. FLAG® Positive Control Lysate is suggested for use as a control in western blotting assays in conjunction with Anti-DYKDDDDK Monoclonal Antibody HRP Conjugate.

IV. Affinity Gel Regeneration

1. Anti-DYKDDDDK Affinity Gel can be regenerated immediately after use by washing with 3 gel volumes of 0.1 M glycine HCl, pH 3.5.
2. The affinity gel must be immediately re-equilibrated with 10 gel volumes of 0.02M Potassium Phosphate, 0.15M Sodium Chloride pH 7.2 with 0.01% (w/v) sodium azide (added as a preservative).
3. Store the used affinity gel at 4°C.

Rockland products are for research use only and are not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Immunochemicals, Inc., P.O. BOX 5199, Limerick, Pennsylvania 19468, USA.

© 2025 Rockland Immunochemicals, Inc. All rights reserved.

FLAG® is a registered trademark of Merck KGaA.