



AccuSignal™  
**Nuclease ELISA Kit**  
**(KJE-4001)**

**Contents:**

- I. Kit Components
- II. Reagents Required but Not Provided
- III. Assay Reagents Preparation
- IV. Sample Preparation & Matrix Effects
- V. Recommended Assay Procedure
- VI. Calculations
- VII. Assay Performance Characteristics
- VIII. Troubleshooting

## Introduction

AccuSignal™ Nuclease ELISA Kit (KJE-4001) is designed to sensitively and robustly quantitate nuclease within therapeutics. This kit is a 96-well microtiter strip format sandwich ELISA that uses pre-immobilized anti-nuclease antibodies, biotin-conjugated detection combined with streptavidin-HRP, and 3,3',5,5'-TMB substrate to measure nuclease in a variety of therapeutic products and materials.

Captured proteins are detected using a biotin conjugated Nuclease antibody combined with HRP conjugated Streptavidin. The addition of TMB develops in a chromogenic change that can be measured photometrically at 450 nm with reference wavelength at 630-650 nm. The OD in the wells is directly proportional to the nuclease concentration of the solution and can be interpolated to calculate the nuclease present in the sample. The nuclease concentration can be calculated based on the 4PL calibration curve generated with the provided nuclease standard.

## I. Kit Components

Table 1. List of Accusignal Nuclease ELISA Kit Components

Component	Item No.	Size
Nuclease Antibody-coated 96-well Strip Plate <sup>1</sup>	KJE-4001B	1 plate
Biotinylated Anti-Nuclease Detection Antibody (100X) <sup>2</sup>	KJE-4601A	110 µL
Nuclease Protein Standard <sup>3</sup>	KJE-0001C	0.2 µg
Streptavidin-HRP (100X) <sup>2</sup>	KJX-0001K	110 µL
Sample Buffer <sup>4</sup>	KJX-0001D	50 mL
Stop Buffer <sup>4</sup>	KJX-0001G	20 mL
TMB Buffer <sup>4</sup>	KJX-0001F	20 mL
Wash Buffer (10X) <sup>2</sup>	KJX-0001E	60 mL
Plate Sealer <sup>4</sup>	KJX-0001H	1 sheet

**Note:** Store kit components at 2°C to 8°C upon arrival. Aliquot the reconstituted protein standard and store at -20 °C.

<sup>1</sup> Antibody-coated 96-well strip plate in a 1 x 8 strip format. Each strip is packaged within a frame enclosed in a resealable foil pouch that includes a desiccant. Ready for use as supplied. Do not use wells if the foil pouch shows evidence of damage, such as tears or holes in the foil pouch.

<sup>2</sup> Dilute prior to use.

<sup>3</sup> Reconstitute and dilute prior to use.

<sup>4</sup> Ready to use as supplied.

## II. Reagents & Materials Required but Not Provided

- Microplate shaker (up to 450 rpm)
- Interval timer
- Multichannel pipettor (50-300 µL)
- Precision single pipettes (10 µL, 35 µL, 100 µL, 1000 µL, etc.)
- Disposable pipette tips
- Deionized water

- Disposable microcentrifuge Tube(s) or microplate
- Polypropylene centrifuge tubes (15 mL)
- Spectrophotometer microplate reader (450 nm absorbance, 630–650 nm reference filter)
- Disposable gloves
- Graduated cylinder
- Reagent reservoirs
- Vortex mixer
- Stir plate & magnetic stir bar
- Absorbent paper

### III. Assay Reagents Preparation

#### Preparation of Standards & Test Samples

1. Reconstitute the standard vial with 1.0 mL deionized water to obtain a final concentration of 200 ng/mL.

**Note:** This is the “reference stock solution” that will be used below to make the standards (see Table 2).

2. Prepare dilutions of standard as indicated in Table 2.
3. Test samples should be diluted in Nuclease Kit Sample Buffer based on empirically determined criteria for each sample.

**Table 2.** Standards, Concentrations, and volumes

Assay Standard	Final Conc. (ng/mL)	Initial Volume (μL)	Add Vol. to Dilution Tubes
A	20.0	900	100 μL of Ref stock solution (RSS at 200 ng/mL)
B	6.9	570	300 μL of A
C	2.4	570	300 μL of B
D	0.8	570	300 μL of C
E	0.3	570	300 μL of D
F	0.1	570	300 μL of E
G	0.03	570	300 μL of F
H	0	1,000	0

## Detection antibody working solution preparation

**Table 3.** Nuclease Detection Antibody Concentration and Volume for a full 96-well assay (1:100 ratio of antibody to buffer)

Nuclease Sample Buffer Volume	11 mL
Biotinylated Anti-Nuclease Detection Antibody Volume	110 $\mu$ L

1. Add 110  $\mu$ L of conjugated antibody to 11 mL of Sample Buffer for use in a full 96-well assay.
2. Mix well by pipette or inversion. Do not vortex.
3. Distribute antibody working solution as described in the assay procedure.

**Note:** Volumes may be adjusted so long as final working concentration remains as specified.

## Streptavidin-HRP Working Solution Preparation

**Table 4.** Streptavidin-HRP Concentration and Volume for a full 96-well assay (1:100 ratio of streptavidin-HRP to buffer)

Nuclease Sample Buffer Volume	11 mL
Streptavidin-HRP Volume	110 $\mu$ L

1. Add 110  $\mu$ L of Streptavidin-HRP to 11 mL of Sample Buffer, respectively for use in a full 96-well assay.
2. Mix well by pipette or inversion. Do not vortex.
3. Distribute Streptavidin-HRP working solution as described in the assay procedure.

**Note:** Volumes may be adjusted so long as final working concentration remains as specified

## Wash Buffer (1X) Preparation

1. Add 50 mL of Nuclease Kit Wash Buffer (10X) to 450 mL of deionized water.

2. Mix for at least 10-minutes using a magnetic stir bar.

**Note:** The wash buffer (1X) can be stored at room temperature (15°C to 25°C) for up to 2-weeks, after which it should be discarded.

## VI. Sample Preparation & Matrix Effect

The performance of the kit was tested in various sample buffers (matrices). Some buffer systems can greatly affect the assay, causing low recovery because of suboptimal conditions (i.e low pH).

### Sample Dilution

Optimal dilution factors for each test sample must be determined by the user. It is recommended to dilute the test sample first to a minimum ratio of 1:2, or ratios as determined by the user, using Nuclease Kit Sample Buffer and perform a serial dilution of the sample with Nuclease Kit Sample Buffer. For low pH buffers, a dilution of 1:10 or up to 1:20 might be required.

**Note:** Remove insoluble protein aggregates by centrifuging the sample.

### Matrix Effect

Validate specific matrix effects if the test sample is not diluted into Nuclease Kit Sample Buffer. Matrix effects can be quantified by spiking known concentrations of Nuclease Protein Standard into test matrix buffers. Matrix effects should be tested with different proportions of the Nuclease Kit Sample Buffer for optimal results. Recovery should be 80-120%

## V. Recommended Assay Procedure

The assay procedure takes approximately 3 hours with a detection range of 0.03–20 ng/mL.

For a reproducible assay we recommend to:

- Mix the samples and reagents before use.
- Perform the washing procedure with the recommended wash buffer volumes.
- Measure all standards and samples in duplicate or triplicate.
- Avoid the use of an automated plate washer unless validated as fit for purpose by the end user.

1. To each well add 100  $\mu\text{L}$  of unknown or standard sample per well and incubate at room temperature for 60-minutes with shaking at 450 revolutions per minutes (rpm) on a shaker.
2. Wash the wells with Wash Buffer as follows:
  - i. Decant the contents of the wells manually with a hard, rapid downward motion. Fluid should be captured in a receptacle designed to collect the waste. Remove all residual reagent from the microplate by tapping it on absorbent paper with the opening facing downwards.
  - ii. Fill each well with 300  $\mu\text{L}$  of Washing Buffer with a multichannel pipettor.
  - iii. Decant the Washing Buffer from the wells with a hard, rapid downward motion. Remove all residual solution from the microplate by tapping it on absorbent paper with the opening facing downwards.
  - iv. Repeat steps ii and iii two more times (total of 3 washings). Do not leave any residual moisture in the wells on each washing step.
3. To each well add 100  $\mu\text{L}$  of detection antibody working solution (see Table 3)

4. Incubate at room temperature for 60-minutes, covered to protect from light, with shaking at 450 rpm.
5. Following 60-minute incubation, wash with Wash Buffer as follows:
  - i. Decant the contents of the wells manually with a hard, rapid downward motion. Fluid should be captured in a receptacle designed to collect the waste. Remove all residual reagent from the microplate by tapping it on absorbent paper with the opening facing downwards.
  - ii. Fill each well with 300  $\mu$ L of Washing Buffer with a multichannel pipettor.
  - iii. Decant the Washing Buffer from the wells with a hard, rapid downward motion. Remove all residual solution from the microplate by tapping it on absorbent paper with the opening facing downwards.
  - iv. Repeat steps ii and iii two more times (total of 3 washings). Do not leave any residual moisture in the wells on each washing step.
6. To each well add 100  $\mu$ L of Streptavidin-HRP working solution (see Table 3)
7. Incubate at room temperature for 20-minutes, covered to protect from light, with shaking at 450 rpm.
8. Following 20-minute incubation, wash with Wash Buffer as follows:
  - i. Decant the contents of the wells manually with a hard, rapid downward motion. Fluid should be captured in a receptacle designed to collect the waste. Remove all residual reagent from the microplate by tapping it on absorbent paper with the opening facing downwards.
  - ii. Fill each well with 300  $\mu$ L of Washing Buffer with a multichannel pipettor.
  - iii. Decant the Washing Buffer from the wells with a hard, rapid downward motion. Remove all residual solution from the microplate by tapping it on absorbent paper with the



- opening facing downwards
- iv. Repeat steps ii and iii two more times (total of 3 washings).  
Do not leave any residual moisture in the wells on each washing step.
  9. Next add 100  $\mu\text{L}$  per well of room temperature TMB solution and incubate the plate with TMB solution at room temperature for 20 minutes (covered and protected from light).
  10. Then add 100  $\mu\text{L}$  of stop solution per well. Gently tap the plate to mix, ensuring no bubbles are formed, and read plate within 5 minutes after stopping the reaction.
  11. On plate reader, measure absorbance at 450 nm with the reference wavelength set at 630–650 nm.

## VI. Calculations

Follow the steps below to estimate the nuclease concentration of the test samples.

1. Calculate the relative OD 450 using the following formula:  
Relative OD 450 = (OD 450 of well) - (OD 630–650 nm of the well)
2. Calculate the mean relative OD 450 of the replicates for each standard solution.
3. Plot the standard solutions data as mean relative OD 450 for each standard solution (Y) vs the respective concentration of the standard solutions (X).
4. Fit the standard solution data with a 4-parameter logistic (4-PL) curve. Weight by  $1/Y^2$  is intended to be used during generation of 4-PL curve
5. Estimate the Nuclease concentration of each test sample well using interpolation from the standard curve. Calculate the average of each respective sample solution concentration.

**Note:** If the assay samples are from dilutions, multiply the

concentrations obtained from interpolations by the dilution factor.

**Note:** If the spectrometer used for the assay does not automatically subtract the reference wavelength, do this manually.

## 9. Assay performance characteristics

Table 5. General kit specifications

Specification	Parameter
Goodness of fit of standard curve (4PL)	$R^2 > 0.98$
Range of the standard curve	1.56 – 200 ng/mL
Precision, intra- and inter-assay CV%	<20%
Sensitivity, lower limit of quantification	3.0 ng/mL
Sensitivity, lower limit of detection	<1.0 ng/mL

## 10. Troubleshooting

Table 6. Troubleshooting

Issue	Possible Solutions
Low Reproducibility	<ul style="list-style-type: none"><li>• Never mix components between kit lots. Each assay is lot-specific and designed for performance when used before the assigned expiration date.</li><li>• Review the washing step procedure. Add recommended volumes of wash buffer to all wells</li><li>• Review preparation and mixing of the samples procedures. Make sure that the samples are well-mixed and not aggregated.</li><li>• Make sure to have proper reagent temperatures throughout the assay.</li><li>• Liquid on the underside of the plate. Carefully dry the underside of the plate with an absorbent towel.</li></ul>

**Table 6.** Troubleshooting

Issue	Possible Solutions
High Background	<ul style="list-style-type: none"><li>• Review washing step procedure. Add recommended volumes of wash buffer to all wells.</li><li>• Review the assay incubation procedure. Follow recommended incubation times.</li><li>• Review preparation and mixing of samples procedure. Make sure that samples are well mixed and not aggregated.</li></ul>
Poor Reactivity or Intra Assay Precision	<ul style="list-style-type: none"><li>• Review preparation and the mixing of samples procedure. Make sure that the samples are well mixed and not aggregated.</li><li>• Review matrix effect: high concentrations of detergents, salts, pH and other matrices characteristics can influence the assay.</li><li>• Use the provided sample buffer to dilute the sample to the minimum recommendations</li></ul>
Poor Fit of the Standard Curve	<ul style="list-style-type: none"><li>• Review washing step procedure. Add recommended volumes of wash buffer to all wells.</li><li>• Review preparation and mixing of samples procedure. Make sure that the samples are well mixed and with no protein aggregates.</li><li>• Review pipetting procedure. Make sure that the samples are pipetted into the bottom of each well. Take care not to touch the bottom of the well with the pipette tip.</li></ul>



Rockland Immunochemicals, Inc.

PO Box 5199, Limerick, PA 19468

+1 484.791.3823

[www.rockland.com](http://www.rockland.com)

Rockland products are for research use only and are not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Immunochemicals, Inc., P.O. BOX 5199, Limerick, Pennsylvania 19468, USA.

© 2023 Rockland Immunochemicals, Inc. All rights reserved.